Instructions for use

DAWID MINI TRANSPORT LIFTING DEVICE
DRVF1
EDITION EJ 01.2022

VITEA[®]
CARE

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ATTENTION

Always follow the basic precautions according to the information included below and regarding safety measures when using the device. Read the entire instructions carefully before use. Keep this user manual for future reference.

1. Introduction

Thank you for choosing our product. We are deeply convinced that it will meet your expectations. Please read this document before first use. This manual was made in order to provide the necessary information on the safe use of the device. The correct application of the recommendations will contribute significantly to extended durability and visual design of the product. We would like to highlight that to date we are constantly improving the quality of our products, which may result in minor changes that are not included in this manual.

2. Product description


Fig. 1. Construction scheme of Dawid mini transport lifting device. The letters are marked with the following structure elements: **A** - base, **B** - frame, **C** - electric actuator, **D** - width adjusting pedal, **E** - control box with battery, **F** - remote control, **G** - lifting arm, **H** - adjusting screw, **I** - hanger, **J** - front castor without brake, **K** - rear castor with brake, **L** - base.

3. Indications

Dawid mini transport lifting device is designed to carry only one patient at a time. For safety during lifting and transport, patient should be under supervision of trained caregiver. Each time in order to maneuverer should use handle. Thanks to the solutions implemented from the device can be used patients affected:

- Patients with paraplegia - including the lower limbs.
- Patients with paralysis.
- Patients with limb loss.
- Patients with cardiovascular failure.
- Patients with reduced joint mobility and muscle weakness.
- Elderly people requiring specialized home or professional care.

4. Contraindications

The mini transport lifting device Dawid equipment is not a device used to transport patients over long distances. Choose the sling recommended by the doctor for the patient. When moving the patient while he is suspended from a carrier attached to the lift, do not move the base wheels through items such as raised carpet frames, door sills, or any uneven surfaces that may cause the lift to shake. If any undesirable effects occur, the patient's transport should be stopped immediately.

5. Usage

- After unpacking check that it contains all the parts in accordance with the list of items above (Fig. 1).
- Make sure that all components do not have visible mechanical damage. In case of doubt, contact the manufacturer.
- Do not use the device until all components are properly connected.
- Maximum permissible user weight - 150 kg
- We are not responsible for any damages resulting from the use of a lift with a higher weight than the permissible one.
- Please do not place the hoist on an inclined surface as it may fall down.
- Do not attempt to exceed any steps, curbs, or other obstacles that could cause injury or injury when the product falls over.
- Do not attempt to repair the battery yourself

5.1. Assembly guide

Step 1. Install the column

Fig. 2. Scheme of mounting the lifting column.

Remove the device and spare parts from the package. Loosen the jam nut and bolt, then raise the frame and align the holes in parallel. Tighten the screws and tighten the nut properly (make sure there is no movement at the frame).

ATTENTION: When mounting the lifting column, avoid standing under it, failure to do so may result in accidental fall of the column and risk of injury.

Step 2. Hoist Arm and Powered Actuator


Fig. 3. Installation of lift arm

Align the bracket with the mast mounting hole, screw in the bolt and washer, install the bottom of the actuator and column in place, and finally install the electric actuator and bracket, tighten all nuts.

Step 3. Install the Controller


Fig. 4. Controller installation scheme

Releasing: Press the release button and pull the battery out in the direction. Installation: Place the battery in place and press the switch in the direction shown in the scheme (Fig. 4). Ensure that the battery cannot be separated without pressing the switch.

Fig. 5. Installation scheme of hooks for the power cord.

Insert the hooks into the side holes of the battery box (Fig. 5) to facilitate storage of the power cord. One end of the power cord should be connected to the control box, while the other end should be connected to a wall outlet to charge the battery. To disassemble the Dawid transport lifting device, proceed in the opposite way than for assembly.

5.2. Operation guide

Electrical controller (remote control):

a) Emergency Stop Button (Emergency Stop)

The emergency stop button interrupts and turns the power off, causing the engine to stop immediately in an emergency.


Fig. 6. Emergency Stop Button operation scheme, positions A and B.

Pressing the Emergency Stop button (red button) (Fig. 6A) will stop the engine. To restore the engine to normal operation, turn the button (Fig. 6B) clockwise to unlock the device. Lifting and lowering the boom are simple and is done using the above pilot / controller. Suspend the remote control after using or adjusting the boom height in a safe and easily accessible place near the driving handle.

NOTE: When the device is not used for more than 3 days, press the Emergency Stop Button to maintain longer battery life.

b) Electrical failure - 'up' and 'down' buttons:

Fig. 7. Schematic of the position of the "up" and "down" buttons.

Use the 'up' or 'down' arrow keys (Fig. 7) to bring the motor to its default position if the controller fails or there is no controller handle.

c) Controller support:

Use the "up" and "down" buttons on the controller handle to control the lift and lowering the lifting arm, after releasing the buttons it will stop moving. If the remote control does not work properly, press the red Emergency Stop Button.

d) LCD display:


Fig. 8. Scheme of controller information displayed, designations: **A** – Over duty on main lifting actuator, **B** – Overload on actuator/s, **C** – Remind to replace main actuator, **D** – Battery level bars, **E** – Error code display.

The display can show tips and information on transport lifting device malfunctions. Battery charge indicator - from the left: 100%, 75%, 50% and 25-0%.

e) Control device connection:


Fig. 9. Electric box controller port diagram.

The device has 4 ports (Fig. 9), from left to right: Handset, lifting actuator (M1), Spread actuator (Optional), Charging plug.

Installation: Insert 4 cables into four different ports one after the other, respectively.

Base width adjustment


Fig. 10. Scheme width adjustment of the base.

The width of the lift base is manually adjusted with a foot pedal. Press the left pedal to extend the base. Press the right pedal to make the base narrower. It is possible to manually lock the base width adjustment by turning the red sleeve (Fig. 10) counter clockwise.

Warning! Locking the wheel during hoist movement will increase the risk of falling down the hoist.

Rear wheel opening and braking


Fig. 11. Foot brake on the rear wheels of the lifter

The rear wheels can be locked to prevent movement. Lock or unlock the rear wheels with your feet. Use your foot to control the brake position. Press the longitudinal plastic element with the foot to block the wheel (Fig. 10A), then press the shorter plastic element to resume movement (Fig. 10B).

WARNING! Locking the rear wheels of the jack increases the risk of the machine tipping over when lifting and lowering. Always keep the support perpendicular to the patient's center of gravity when working.

Charging (power supply)

When the battery has only one line left, the control box emits 4 beeps and remind you to immediately charge. The engine will not run without charging.

Battery charging procedure

When you hear a beep, turn off the emergency stop button and connect the control box. When the battery is fully charged, the battery charge status shows 4 dashes (100%) and the LED on the control panel does not flash.

WARNINGS:

1. Turn off the emergency stop button and do not use the elevator while charging.
2. Take out the control box and keep it in a shaded place when the lift will not be used for a long time. Charge the battery once a month.
3. Do not charge the battery in a humid atmosphere.
4. In the event that the lift is not in daily use, it is recommended to disable the emergency stop button to cut off the current flow and maintain the battery charge.

Note: During charging, the emergency stop button must be turned off and charge, the hoist will be switched off during charging. When the jack is not used for a long time, it must be

reused to maintain the battery's Active Lock performance. Do not charge the lift in high humidity.

6. Safety of use

- a) Before using the lifting device refer to this manual and keep it for future reference.
- b) The lifting device is designed to use indoors only on hard and even surfaces.
- c) Using the lifting device on an unstable surface (covered with ice, wet, greasy, sandy) is unacceptable- it causes risk of falling,
- d) You can only use the lifting device on surfaces where all four wheels are touching the ground.
- e) Do not leave the lifting device with the user close to open sources of fire.
- f) Be aware that the surface temperature of the frame parts may be increased after exposure to the sun or other heat source. Furthermore, with low temperatures there is a danger of overcooling metal parts of lifting device. It may cause burns or frostbite parts of the body. g) In order to transfer lifting device, grab the frame. Never transfer the lifting device by grabbing movable elements or control unit.
- h) **The lifting device is designed for people up to 150 kg.**
- i) Lifting device must be easy to fold and to unfold. In case of problems with folding, you should re-read the instruction or contact the manufacturer.
- j) Brakes cannot slow down the lifting device during ride, you should use them only to prevent the unintended movements.
- k) Oil, grease, water or other chemicals, may cause that brakes may have reduced effectiveness.
- l) It is extremely important to keep all lifting devices components clean in order to maintain the product parameters.
- m) Storing the lifting device in areas with very high humidity (baths, laundry, etc.) or in the open air, also the lack of hygiene causes rapid loss of operating parameters and aesthetic of the product.
- n) It is unacceptable to use the lifting device after discovering any visible or detectable damages or when the wear of some parts may cause an accident. All the repairs need to be performed by skilled personnel and during the warranty period only by the manufacturer's service.
- o) During transport (e.g. In the trunk), do not burden lifting device elements, it may weaken or damage its structure.
- p) Keep out of the reach of children under five years of age. q) Should not insert fingers between the lifting device components - risk of injury or entrapment.
- r) Keep extreme caution while unfolding / folding the lifting device. Lack of attention may result in finger entrapment.
- s) Always maneuvers the lifting device with the handle provided.
- t) Do not use a sling that is not recommended for the lift.
- u) Only trained and qualified caregivers should transfer a person.

- v) Lifting device is designed for transferring patients only. Do not use the lifting device for any other purpose.
- w) Before each patient lifting, make sure that base legs were blocked in maximum wide position.
- x) While moving the patient, should bypass obstacles at safe distance, which might cause the lifting device to tip over.
- y) Before each lifting device usage, make sure that slings are in good condition. Usage slings showing signs of damage is forbidden.
- z) Before using lifting device, should charge the battery for about 6 – 8 hours. Make sure, the red warning lights is off.

7. Everyday use and maintenance

7.1. Lifting the patient

- First, set the ride base in the maximum enlarged position. To do this, press the pedal of the head adjusting the width of the transport lifting device.
- For easy and comfortable lifting of the patient should be prepared a transporter lifting device, a sling with a leg separator, a sling with head stabilization and a wheelchair.
- The patient should lie in the middle of the bed, then turn he/she to the side, in the opposite direction to the caregiver.
- The sling should be placed under the patient.
- Examples of slings are shown in the figures 12 & 13.
- Turn the patient towards the caregiver and place he/she in the middle of the sling so that the back loops are in the upper back of the patient and the thigh loops be placed within the patient's thighs.
- Hanging loops should be hung on a four-point cradle.
- Loops have different colors to make it easier to put on slings at the same height.

Fig. 12. Sling with a leg separator

Fig. 13. Sling with head stabilization

- Make sure hanging loops are tight and make sure the patient is in a comfortable position.
- Lock the rear wheel brakes and lift the patient by pressing the corresponding button on the remote control.
- Should lift the patient until his feet are not touching the bed or the floor.
- Then unlock the rear wheel brake and carry the patient over the wheelchair.
- The brakes should be locked both in the transport lifting device and in the wheelchair.
- The next step is to press the "down" button on the remote control to gradually lower the patient.
- When descending, help the patient to take a correct and a comfortable sitting position.
- After setting down the patient in a wheelchair, remove the sling loops from the four-point cradle.
- If using a sling with a leg separator, make sure that the patient is properly dressed before moving the lift towards the movable toilet chair.
- Make sure that the hanging loops are tight and that the patient is in a comfortable position.

The patient can be placed on the sling in a sitting position, e.g. when transported from a wheelchair or bed as follows.

- Place the sling around the back of the patient.
- Donning is easier if the patient bends slightly forward.
- Guide the sling carefully at the patient up on the height of the tailbone.
- Move the thigh belts forward at patient's thighs until they are on both sides at the same height.
- Be careful not to bend the material.

- Wear the seat belts under the patient's slightly parted legs.
- Donning is easier if the patient slightly raises thighs.

In order for the patient to sit comfortably on a wheelchair, it should be planted as deep as possible.

- During the process, lowering the patient, lightly press the patient's knee so that his/her back is closest to the backrest of the wheelchair.

Slings:

Clean the seats with a damp cloth. Difficult dirt should be removed with a commercially available mild cleaner. Stains can be removed with a cloth or a soft brush. Wash with a wet cloth (with clean water) and allow to dry. Before washing, follow the instructions on the label of the sling.

<p>Sling cleaning instructions:</p> <p>Do not use:</p> <ul style="list-style-type: none"> • Solvents. • Toilet cleaners. • Sharp brushes and hard objects. • Caustic cleaners. 	<p>Washing instructions:</p> <p>The belts can be machined washed at 30 ° C.</p> <ul style="list-style-type: none"> • Do not use chlorine based bleaching agents. • Do not tumble dry. • Do not iron
---	---

Foldable and portable frame:

- Lower the control unit and then fully attach the lift frame.
- Remove the screw before folding the frame.
- Release and support the plate before raising the frame
- The scheme of folding the frame is shown in figure 14.


Fig. 14. Scheme of folding the elevator frame.

7.2. Maintenance and conservation

Routine inspection listed below shall be performed every day for the sake of safety.

1. Check if there is any damage or wear by visual inspection.
2. Make sure screws and nuts on lifting arm are secured.
3. Check if hook on lifting arm works properly.
4. Controller, actuator and base width adjustment verification.
5. Emergency descend function verification (both manually and electrically)

If needed, to clean the patient lift with warm water or alcohol. Make sure there is no dust and hair on the wheels.

Attention: to clean the patient lift with chemicals with Phenol and chlorine is prohibited, which may cause damage on plastic powder or plastic material.

Transportation and Storage:

Press the emergency stop button for transportation or not using the patient life for a certain of time. Temperature for storage of patient life shall be higher than 0°C and relative moisture ≤ 60%.

Inspection:

Make a full inspection at least once a year. Pay close attention on those worn part such as connection and screw bolt.

7.3. Troubleshooting




Issue description	Troubleshooting
<p style="text-align: center;">Patient lift cannot work (ascend and descend)</p>	<ol style="list-style-type: none"> 1. Check if the emergency stop button is pressed down. 2. Verify the controller and cables from actuator is properly connected. 3. Check the battery level from display screen on the control box. 4. Verify if there is any problem or damage on battery contact panel. 5. Ask for after-sale service if patient lift still cannot work.

Charger cannot charge the battery	<ol style="list-style-type: none"> 1. Check if the charger fully plugs into socket. 2. Verify if there is any problem or damage on battery contact panel. 3. Ask for after-sale service if patient lift still cannot work.
Patient lift stops when ascending	<ol style="list-style-type: none"> 1. Check if the emergency stop button is pressed down. 2. Employ the emergency descend device. 3. Ask for after-sale service if patient lift still cannot work.

8. Technical specifications

Technical parameter	Parameter value
Frame	Steel, powder coating
Chassis length	950 mm
Chassis width	520-580 mm
Lifting height	1230 mm
Handle height	780-1870 mm
Castor height	115 mm
Motor	DC 24V 8000N
Battery	DC 24V 5AH, interchangeable battery Approx. 60 lifts per full charge
Charger	24V-1A
Charging time	5 hours per time
Front castor	3" solid tires
Rear castor	4" solid tires with locking
Push handle	Push handle with rounded corner
Max. user weight	150 kg
Addons	Standard sling or Bathing sling
	Standard Mode / Stand Aid Mode Interchangeable
Package	1300x570x350 mm
G. / N. W.	46/40 KGS

9. Symbols

	Confirmation of product compliance with EU standards
	Production date
	Manufacturer
	See the instructions for use
	WARNING
	Catalogue number
	Serial number
	Protect from moisture
	Protect from sunlight
	Permissible user weight
	Indoor use only
	Possibility of fingers jamming
	Temperature of storage and use
	Medical device

DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT

This symbol on the products or accompanying documents means that used electrical or electronic products should not be disposed of with the normal household waste. For proper disposal, refurbishment or recycling, please take these products to the collection points for this type of waste. Alternatively, in some EU countries or other European countries, you may return your product to your local retailer when purchasing a similar new product. By disposing of this product correctly, you will help to conserve valuable natural resources and support the prevention of potential negative effects on the environment and human health, which could result from incorrect disposal.

For further information, please contact your local municipality or collection point waste collection points. Penalties may be imposed in accordance with local regulations if this product is disposed of in an inappropriate manner. For parties within the European Union, If you wish to discard electrical or electronic equipment, please obtain the necessary information from your retailer or supplier.

Elimination in countries outside the European Union. This symbol is valid in the European Union. If you wish to decommission this appliance, please obtain the correct decommissioning information from your local authorities or from your retailer.

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