


2991041

V200 DRIVE SYSTEM HEM2				Plan No. E000651 / F
Sheet 1/2				
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	21/03/2019	21/03/2019	
Initials	DD	SD	PJ	

E000651 / F

Material			
No.	Plan	Article	Quantity
1	X002976	1995343	2
2A	Y021225	2012800	1
2B	Y022137	2025100	1
3	Y017173	1742454 (R)	1
4	Y017205	1742554 (L)	1
5	M6x40 ISO7380	1006111	4
6	M6 DIN986	1008492	4
7	---	1901208	1
8	---	1901745	1
9	---	1994568	4
10	DIN471	1995259	4
11A	Y021227	2012900	1
11B	Y021231	2013000	1
11C	Y022136	2025200	1
12A	X005892	G103100	1
12B	X005891	G103000	1
13	X005893	1915848	1
15	M6 DIN7349	1003321	1
16	M6x35	2990215	1

2A Seat width 38 to 52

2B Seat widths 53 to 60

11A Seat widths 38 to 44

11B Seat widths 46 to 52


11C Seat widths 53 to 60

12A Seat widths 38 to 44

12B Seat widths 46 to 60

2991041

E000651 / F

V200 DRIVE SYSTEM HEM2				Plan No.
Sheet 2/2				E000651 / F
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	21/03/2019	21/03/2019	
Initials	DD	SD	PJ	